5

10

15

20

25

What is Claimed is

- 1. (amended) A method for measuring the amount of β ig-h3 protein comprises the following steps:
 - 1) Preparing recombinant β ig-h3 proteins comprising 4th fas-1 domains, their fragments or derivatives, as standard proteins;
 - 2) Preparing specific ligands against the above recombinant proteins, their fragments or derivatives of the above step 1; and
 - 3) Measuring the amount of β ig-h3 protein of samples with the method using binding reaction of ligands of the above step 2 with the recombinant proteins, their fragments or derivatives of the above step 1.
- 2. The method for measuring the amount of β ig-h3 protein as set forth in claim 1, wherein the ligands of step 1) are selected from a group consisting of antibodies, RNA, DNA, lipids, proteins, organic compounds and inorganic compounds.
- 3. The method for measuring the amount of β ig-h3 protein as set forth in claim 1, wherein the

PCT/KR2002/001975 IPEA/KR 09. 09. 2004

specific binding reaction of step 3) is antigenantibody reaction.

4. The method for measuring the amount of β ig-h3 protein as set forth in claim 3, wherein the antigen-antibody reaction is performed by a. selected from method a· group consisting immunoblotting, immunoprecipitation, ELISA, RIA, protein chip, rapid assay and microarray.

10

5

5. (amended) The method for measuring the amount of β ig-h3 protein as set forth in claim 3, wherein antigen-antibody reaction of step 3) comprises the following steps:

15

1) Coating recombinant β ig-h3 proteins comprising fas-1 domains, their fragments derivatives to matrix;

20 -

2) Reacting antibody against the protein of the above step 1, its fragments or derivatives with sample;

3) Adding the reactant of the above step 2 to the coated protein of step 1 and waiting for reaction, and then washing thereof; and

25

4) Adding the secondary antibody to the reactant of the above step 3 for further reaction, and then measuring OD.

5

20

25

- 6. The method for measuring the amount of β ig-h3 protein as set forth in anyone of claim 1-5, wherein the β ig-h3 protein is human β ig-h3 protein having amino acid sequence represented by SEQ. ID. NO 3 or mouse β ig-h3 protein having amino acid sequence represented by SEQ. ID. No 5.
- 7. (amended) The method for measuring the amount of β ig-h3 protein as set forth in anyone of claim 1-5, wherein the recombinant β ig-h3 proteins comprising 4^{th} fas-1 domains have 1-10 repeatedly-linked fas-1 domains.
- 8. The method for measuring the amount of β ig-h3 protein as set forth in claim 7, wherein the fas- 1 domain of β ig-h3 is selected from a group consisting of sequences represented by SEQ. ID. No 7, No 8, No 9 and No 10.
 - 9. The method for measuring the amount of β ig-h3 protein as set forth in claim 1, wherein the sample can be any body fluid including urine, blood or synovial fluid.
 - 10. A diagnostic kit for the renal diseases, hepatic

PCT/KR2002/001975 IPEA/KR 09. 09. 2004

diseases, rheumatoid arthritis or cardiovascular diseases comprising β ig-h3 protein or recombinant proteins of fas-1 domain in the β ig-h3 protein (including their fragments or their derivatives) and their ligands.

Car South

5

10

20

25

- 11. The diagnostic kit as set forth in claim 10, wherein the ligand is selected from a group consisting of antibody specifically binding to β ig-h3 protein, fas-1 domain of β ig-h3, their fragments or derivatives, RNA, DNA, lipids, proteins, organic compounds and inorganic compounds.
- 15 12. The diagnostic kit as set forth in claim 11, wherein the ligand is antibody.
 - 13. The diagnostic kit as set forth in claim 12, wherein the kit additionally includes buffer solution, secondary antibody, washing solution, stop solution or coloring substrate.
 - 14. The diagnostic kit as set forth in claim 10, wherein the β ig-h3 protein is human β ig-h3 protein having amino acid sequence represented by SEQ. ID. No 3 or mouse β ig-h3 protein having

amino acid sequence represented by SEQ. ID. No 5.

15. The diagnostic kit as set forth in claim 10, wherein 1 or 2-10 4^{th} fas-1 domains of $\beta\,\text{ig-h3}$ protein are repeatedly linked.

5

10

16. The diagnostic kit as set forth in claim 15, wherein the fas-1 domain of β ig-h3 is selected from a group consisting of sequences represented by SEQ. ID. No 7, No 8, No 9 and No 10.